

Prosthetic armamentaria and auxiliary components

Hex Tools

1.25mmD Hex Tools for Abutment Screws and Fixation Screws



HX1.25 HXL1.25 THX1.25 THXL1.25 0721

Removal Tool

for Internal Hex
3.8mm D Implant Components



TLRT2

Removal Tool

for Internal Octagon
4.8mm D Implant Components



HLRTX2

Torque Wrenches and Inserts



TW30

30 Ncm Torque Wrench used to tighten all components and screws attaching directly into the implant



TW1.25



TW20

20 Ncm Torque Wrench used to tighten screws attaching directly into an abutment.



TW1.25L

Reamers for "Cast-To" or Castable Components




Reamer for Copings PR

Reamer for SPGA MRI

Prosthetic Kit Tray - 2320



Cap Attachment Instruments - CAI



Mandril for Castable Ball Pattern


Nylon Liner Insertion Tool

Nylon Liner Reaming Tool


Cap Attachment System - CAS



Dolder Gold Bar System - DGB



Round Gold Bar System - HGB



Hader Clip Bar System - BS1

